



## Mucosa resection snare

### Model „EVK“ Düsseldorf

Symmetrical flat feather snare

- monofil flat wire for stabilization in the rear area
- twisted wire in the front cutting area





ARTICLE NUMBER SYMMETRICAL FLAT FEATHER SNARE

| Article no.      | Diameter mm | Length cm | Opening width x y in mm | PU |
|------------------|-------------|-----------|-------------------------|----|
| 99 05 20 11 21 6 | 2,3         | 230       | X = 15 / Y = 35         | 5  |
| 99 05 20 11 22 6 | 2,3         | 230       | X = 20 / Y = 40         | 5  |
| 99 05 20 11 23 6 | 2,3         | 230       | X = 30 / Y = 60         | 5  |

We reserve the right to make technical changes and to deviate from the depicted images.

